In case of complications and at the discretion of the surgeon, the implant can be re-used using conventional removal techniques.

WARNINGS AND PRECAUTIONS

Precautions for use

- Risk of interference with medical imaging: MRI/CT scan: the patient should be instructed to systematically declare if he has undergone a surgical procedure.
- The safety and use of the device has not been assessed in a magnetic resonance (MR) environment.
- Heating or migration of the device has not been tested in a magnetic resonance (MR) environment.

Warnings

- The device is not designed to withstand an immediate load after surgery and does not allow for immediate resumption of the patient's activities. If necessary, immobilise during osteosynthesis.
- Device is intended for single use, do not reuse the device to avoid cross contamination.
- Do not re-sterilize the device.

SIDE EFFECTS

- Pain
- Pseudarthrosis or non-union.
- Hypersensitivity to vanadium and/or aluminium.

Please report any serious incident that occurred in relation to the device to Novastep (regulatory@novastep-ortho.com) and to the local competent authority for medical devices.

STORAGE CONDITIONS

Store in a dry place, at ambient temperature.

DISPOSAL

The explanted device is in direct contact with biological tissue and should be disposed of in a biohazardous waste bin according to hospital protocol.

There are no recommendations for the disposal of packaging or any unused device. Refer to the applicable laws for product recycling.

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CONTRA-INDICATIES

- Severe muscular, neurological or vascular deficiency in the extremity concerned.
- Bone destruction or poor bone quality, likely to impair implant stability.
- Hypersensitivity to vanadium and/or aluminium.

DIRECTIONS FOR USE

CAUTION: It is at the discretion of the practitioner to choose the product according to the patient's clinical condition and current surgical practices.

CAUTION: the implant is supplied sterile and must be used in a sterile environment only to avoid contamination risk.

- Select the appropriate device for the indication.
- Before use check:

- the condition of the packaging, if the packaging is damaged, do not use the device and dispose of it according to the disposal section.
- the shelf-life of the device, if the shelf-life is out of date, return the product to Novastep.
- Implant the device using conventional technique.

- CAUTION:*** centrolock® is only compatible with the surgical instruments supplied by Novastep, do not use other instruments from other manufacturers. Refer to the surgical technique (**Air-ST, Air-Ank-CUT, Air-Ank-ST01**) for the list of accessories compatible with the device.

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