Manufactured by:

Stryker GmbH  
Bohnackerweg 1  
2545 Selzach, Switzerland  
Phone: +41 (0) 32 641 66 66  
Fax: +41 (0) 32 641 66 60

CE 0123

V15185 Rev. B, 08-2015 (19-0349)

EN

#### ALL OF THESE INSTRUCTIONS FOR USE MUST BE READ CAREFULLY PRIOR TO CLINICAL USE

Instructions for Use – MCP, MCP-X and PIP (Metacarpophalangeal Prosthesis, Pre-Flex Metacarpophalangeal Prosthesis and Proximal Interphalangeal Prosthesis)

**CAUTION: United States Federal Law restricts this device to sale, distribution and use by or on the order of a physician.**

#### DESCRIPTION

These devices are fabricated from medical grade silicone elastomer. In some implants, this material is internally reinforced with medical grade polyester mesh and/or fibers. Suture ties are provided in some implants for temporary stabilization.

Polyester mesh or velour external surfaces, when present, promote long term stabilization through tissue ingrowth.

#### INDICATIONS

MCP and MCP-X

Intended for Replacement of the Metacarpophalangeal joint of the hand which has been damaged by rheumatoid, osteo or post traumatic arthritis.

#### PIP

Silicone elastomer implants are designed to assist in joint arthroplasty. Clinical indications for use of these devices include:

1. Rheumatoid arthritis
2. Osteoarthritis
3. Ankylosed joints or those with limited range of motion which have not responded to conservative treatment
4. Nonfunctional joint due to inadequate bony alignment and joint space which cannot be restored by soft tissue reconstruction alone
5. Destroyed articular surface(s).

#### CONTRAINDICATIONS

MCP and MCP-X

General contraindications for the use of these implants for joint reconstruction include:

1. A patient with skin, bone, circulatory and/or neurological deficiency
2. Nonfunction and irreparable musculotendinous system
3. Active sepsis
4. Inadequate bone stock.

#### PIP

1. A psychologically unsuitable patient
2. A patient with skin, bone, circulatory and/or neurological deficiency
3. Nonfunctioning and irreparable musculotendinous system
4. Active sepsis
5. Inadequate bone stock.

#### WARNINGS (See also the Patient Counseling Information Section),

Strenuous loading, excessive mobility, and articular instability all may lead to accelerated wear and eventual failure by loosening, fracture, or dislocation of the device. Patients should be made aware of the increased potential for device failure if excessive demands are made upon it.

Implants are mechanical devices that can be worn away, fatigued or broken. An implant site may become infected, painful, swollen, or inflamed. The status of the adjacent bone and soft tissue may be inadequate to support the implant, or may deteriorate in time resulting in instability, deformity, or both. The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time,

necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon. Single use is defined as use of one implant or instrument on a single patient in a single surgical procedure. Reuse of instruments designated as single use has been associated with necrosis of bone leading to implant failure, sepsis and/or communication of potentially lethal viruses.

#### PRECAUTIONS

The implant is provided sterile in an undamaged package. If either the implant or the package appears damaged, expiration date has been exceeded, or if sterility is questioned for any reason, the implant should not be used. **Do not re-sterilize**. Implants should be kept immersed in sterile saline solution after removal from sterile wrapping prior to implantation. Implants should only be handled with blunt instruments to avoid scratching, cutting or nicking the device. Meticulous preparation of the implant site and selection of the proper size implant increase the potential for successful outcome. A complete range of trial sizers are available to aid in bone preparation.

It is suggested that the proper size implant be removed from its sterile package only after the implant site has been prepared and properly sized.

#### ADVERSE EVENTS

Potential adverse events reported with joint prostheses include pain, loosening, fracture, dislocation, or infection. There have been some reports of patients with silicone sensitivity reactions following joint replacement.

Implantation of materials such as silicone may result in foreign body reaction adjacent to the implant site. Injury to the surrounding nerves, blood vessels, tendons, or soft tissues can occur as a consequence of implanting this device. In some patients, wear particles from silicone elastomer implants used in bone and joint reconstruction may participate in, or exacerbate, synovitis or bone cyst complications in contiguous bone. Contributing factors have been reported to include the use of implants in physically overactive patients, associated preoperative pathology such as cysts and degenerative changes, intraoperative temporary stabilization with K-wires, subluxated implants, implant over or undersizing, and uncorrected or recurrent deformity.

#### SURGICAL PROCEDURES

An operative technique for use of this implant device is available. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

#### PATIENT COUNSELING INFORMATION (See also Warnings)

Providing each patient scheduled for implant surgery with documented counseling of potential complications and alternatives, which may include non-implant procedures such as soft tissue reconstruction or arthrodesis, prior to surgery is necessary. In addition to the patient related information contained in the Warnings and Adverse Events sections, the following information should be conveyed to the patient: 1) Adverse effects may necessitate reoperation, revision, or fusion of the involved joint. 2) While the expected life of joint replacement components is difficult to estimate, it is finite. These components are made of foreign materials which are placed within the body for the potential restoration of mobility or reduction of pain.

However, due to the many biological, mechanical and physicochemical factors which affect these devices, the components cannot be expected to withstand the activity level and loads of normal healthy bone for an unlimited period of time.

#### STERILIZATION

This component has been sterilized by ethylene oxide. The implant is provided sterile in an undamaged package. **Do not re-sterilize**. Trial sizer components are available to avoid having to open the sterile package prior to prosthesis implantation. The implant should be removed from its sterile package only after the implant site has been prepared and properly sized.

For symbol explanations, please refer to the table at the end of this document

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