fr	Agrafe d'ostéosynthèse compressive
en	Compressive osteosynthesis staple
de	Osteosyntheseklammer
nl	Compressieklem voor osteosynthese
es	Grapa de compresión para osteosíntesis
it	Cambre di compressione per osteosíntesi
pt	Grampo de osteossintese compressiva
sv	Kompressivt osteosyntes
da	Komprimérende osteosyntese
fi	Puristava osteosynteesi niitti
no	Kompressivt osteosyntes
pl	Zszwka kompresyjna do osteosyntez
el	Συσκευή συμπιεστικής οστεοσύνθεσης

(CE): 2014-06

novastep® S.A.S.
2, Allée Jacques Frimot
35000 Rennes - France
Tel.: +33 (0)2 99 33 86 50
Fax : +33 (0)9 70,29,18,95
contact@novastep-ortho.com

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[en] Osteosynthesis compressive staple intended for the fixation of osteotomy and arthrodesis of the foot.

PRESENTATION

Each (double-wrapped) pack contains a ready-to-use sterile implant with foam holder.

Gamma sterilisation.

COMPOSITION

Nitinol alloy | ASTM F2063

Element Composition (w/w) %

Nickel 54,5 to 57,0

Carbon, max 0,040

Cobalt, max 0,050

Copper, max 0,010

Chromium, max 0,010

Hydrogen, max 0,005

Iron, max 0,050

Niobium, max 0,025

Nitrogen, max 0,005

Oxygen, max 0,040

Titanium Balance

TARGET GROUPS

Target population

Adult patients

Intended user

Orthopaedic surgeons

U.S. - CAUTION: American federal laws require this device to be supplied by a physician or in accordance with a medical practitioner's prescription.

INDICATIONS FOR USE

Osteosynthesis compressive staples are indicated for the correction of bone deformities (Hallux valgus, Lesser ray deformities, Hammertoe deformity, Coalitions, Rheumatoid foot deformity, Cavovarus foot, Flatfoot), and treatment of foot joint stiffness (Hallux rigidus, Foot arthritis).

U.S. only: The compressive osteosynthesis staples are indicated for hand and foot bone fragements osteotomy fixation and joint arthrodesis.

CONTRA-INDICATIONS

• Severe muscular, neurological or vascular deficiency in the extremity concerned.

• Bone destruction or poor bone quality, likely to impair implant stability.

• Hypersensitivity to Nitinol.

DIRECTIONS FOR USE

CAUTION: It is at the discretion of the practitioner to choose the product according to the patient's clinical condition and current surgical practices.

CAUTION: the implant is supplied sterile and must be used in a sterile environment only to avoid contamination risk.

• Select the appropriate device for the indication.

• Before use check:

- the condition of the packaging, if the packaging is damaged, do not use the device and dispose of it according to the disposal section.

- the shelf-life, if the shelf-life is out of date, return the product to Novastep.

• Implant the device using conventional techniques.

CAUTION: arcad® is only compatible with the surgical instruments supplied by Novastep, do not use other instruments from other manufacturers. Refer to the surgical technique (Arc-ST) for the list of accessories compatible with the device.

Removal procedure
In case of complications and at the discretion of the surgeon, the implant can be removed using conventional removal techniques.

WARNINGS AND PRECAUTIONS

Precautions for use

- Risk of interference with medical imaging: MRI/CT scan: the patient should be instructed to systematically declare if he has undergone a surgical procedure.
- The safety and use of the device has not been assessed in a magnetic resonance (MR) environment.
- Heating or migration of the device has not been tested in a magnetic resonance (MR) environment.

Warnings

- The device is not designed to withstand an immediate load after surgery and does not allow for immediate resumption of the patient's activities. If necessary, immobilise during osteosynthesis.
- Device is intended for single use, do not reuse the device to avoid cross contamination.
- Do not re-sterilize the device.

Side effects

- Pseudarthrosis.
- Infection, bruising, venous thrombosis, pulmonary embolism, cardiovascular problems.
- Hypersensitivity to Nitinol.

Please report any serious incident that occurred in relation to the device to Novastep (regulatory@novastep-ortho.com) and to the local competent authority for medical devices.

Storage conditions

Store in a dry place, at ambient temperature.

Disposal

The explanted device is in direct contact with biological tissue and should be disposed of in a biohazardous waste bin according to hospital protocol.

There are no recommendations for the disposal of packaging or an unused device. Refer to the applicable laws for product recycling.

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