Static Staple

IFU0005
Rev J / 2023-12-18

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[en] Osteosynthesis static staples intended for the fixation of mono-cortical phalangeal osteotomy of the foot (Akin osteotomy).

PRESENTATION

Each (double-wrapped) pack contains a ready-to-use sterile implant. Gamma sterilisation.

COMPOSITION

Stainless steel | ASTM F138 | ISO 5832-1

TARGET GROUP(S)

• Target population: Adult patients

• Intended user: Orthopaedic surgeons

INDICATIONS FOR USE

Osteosynthesis static staples indicated for Hallux Valgus.

CONTRA-INDICATIONS

• Severe muscular, neurological or vascular deficiency in the extremity concerned.

• Bone destruction or poor bone quality, likely to impair implant stability.

• Hypersensitivity to nickel and/or chrome.

DIRECTIONS FOR USE

CAUTION: It is at the discretion of the practitioner to choose the product according to the patient's clinical condition and current surgical practices.

CAUTION: the implant is supplied sterile and must be used in a sterile environment only to avoid contamination risk.

• Select the appropriate device for the indication.

• Before use check:

- the condition of the packaging, if the packaging is damaged, do not use the device and dispose of it according to the disposal section.

- the shelf-life of the device, if the shelf-life is out of date, return the product to Novastep.

• Implant the device using conventional techniques.

CAUTION: static staple is only compatible with the surgical instruments supplied by Novastep, not with other instruments from other manufacturers. Refer to the surgical technique (SSS-ST02) for the list of accessories compatible with the device.

Removal procedure

In case of complications and at the discretion of the surgeon, the implant can be removed.

WARNINGS AND PRECAUTIONS

Precautions for use

• Risk of interference with medical imaging: MRI/CT scan: the patient should be instructed to systematically declare if he has undergone a surgical procedure.

• The safety and use of the device has not been assessed in a magnetic resonance (MR) environment.

• Heating or migration of the device has not been tested in a magnetic resonance (MR) environment.

Warnings

• The device is not designed to withstand an immediate load after surgery and does not allow for immediate resumption of the patient's activities. If necessary, immobilise during osteosynthesis.

• The device is intended for single use, do not reuse the device to avoid cross contamination.

• Do not re-sterilize the device.

SIDE EFFECTS

• Pseudarthrosis.

• Infection, bruising, venous thrombosis, pulmonary embolism, cardiovascular problems.

• Hypersensitivity to nickel and/or chrome.

Please report any serious incident that occurred in relation to the device to Novastep (regulatory@novastep-ortho.com) and to the local competent authority for medical devices.

STORAGE CONDITIONS

Store in a dry place, at ambient temperature.

DISPOSAL

The explanted device is in direct contact with biological tissue and should be disposed of in a biohazardous waste bin according to hospital protocol.

There are no recommendations for the disposal of packaging or an unused device. Refer to the applicable laws for product recycling.

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