MANUFACTURED WITH
PEEK ALTERA®
SHAPE MEMORY BIOPOLYMER

ECLIPSE™ SOFT TISSUE ANCHOR

Instructions for Use

English (ENG)

1. DEVICE DESCRIPTION
The Eclipse™ Soft Tissue Anchors, consisting of non-absorbable PEEK Altera® material, are devices for use in the fixation of ligaments, tendons, or soft tissue grafts to bone in shoulder, elbow, knee, foot/ankle, and hand/wrist procedures. The Eclipse Retention Rod and Drive Assembly are manufactured from stainless steel. The optional suture is USP No. 0 Force Fiber Suture.

The Eclipse Sheath, Eclipse Bullet, and the Eclipse Disposable Gun together comprise the Eclipse Device.

2. INDICATIONS FOR USE
The Eclipse Soft Tissue Anchor is intended for fixation of soft tissue to bone in the shoulder, foot/ankle, knee, hand/wrist, and elbow in the following procedures:

- SHOULDER - Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsule Shift or Capsulolabral Reconstruction
- FOOT/ANKLE - Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Midfoot Reconstruction, Metatarsal Ligament Repair, Flexor Hallucis Longus for Achilles Tendon Reconstruction, Tendon Transfer in the Foot and Ankle
- KNEE - Anterior Cruciate Ligament Repair, Posterior Cruciate Ligament Repair, Meniscus/Medial/Lateral Ligament Repair, Lateral Collateral Ligament Repair, Patellofemoral Repair, Quadriceps Tendon Repair, Iliotibial Band Tenodesis
- ELBOW - Biceps Tendon Reattachment, Ulnar or Radial Collateral Ligament Reconstruction
- HAND/WRIST - Scapholunate Ligament Reconstruction, Ulnar Collateral Ligament Reconstruction, Radial Collateral Ligament Reconstruction, Carpalometacarpal Joint Arthroplasty (Basal Thumb Joint Arthroplasty), Carpal Ligament Reconstructions and Repairs, Tendon Transfer in the Hand/Wrist

3. CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS section of this guide.
- Patients with an active local or systemic infection.
- Conditions which tend to result in healing such as blood supply limitations or previous infections.
- Skeletally immature patients where the implanted device would cross open epiphyseal plates.
- Insufficient quality or quantity of bone, comminuted bone surfaces, or pathologic conditions such as cystic change or severe osteopenia that would impair the ability of the Eclipse Soft Tissue Anchor to securely fixate to the bone.
- Inadequate neuromuscular status (e.g. paraparesis, inadequate muscle strength).
- Patients with conditions such as mental illness, senility or alcoholism that tend to restrict his or her ability or willingness to follow postoperative instructions during the healing process.
- Patients with foreign body sensitivity, suspected or documented material allergy or intolerance. Where material sensitivity is suspected, appropriate tests should be conducted and sensitivity ruled out prior to implantation.

4. WARNINGS

- The Eclipse Soft Tissue Anchor is supplied sterile for single use only.
- DO NOT ATTEMPT TO RESTERILIZE. Resterilization may result in loss of proper mechanical function of the device and could result in patient harm.
- Carefully inspect product packaging and all device components for damage or defects prior to use. Do not use any device if it appears defective, damaged or otherwise compromised.
- Do not modify the implant. Modified devices may not perform as intended and could result in patient injury.
- The Eclipse Soft Tissue Anchor is not approved for fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine. Use of the implant or system components for these purposes may result in patient injury.

5. PRECAUTIONS

- The Eclipse Soft Tissue Anchor should be used only by those physicians who have been trained in the appropriate, specialized procedures. Knowledge of appropriate surgical techniques, instrumentation, proper selection and placement of implants, and postoperative patient care and management are essential to a successful outcome.
- Correct selection of the right Eclipse Anchor Components is extremely important. Carefully select the appropriate device size based on the needs of each individual patient.
- Always handle the Eclipse Soft Tissue Anchor carefully. The surface of the implant must always be protected from damage during handling. Avoid contacting the implant with other tools or materials that could notch, scratch or otherwise damage the implant surface. Damage to the implant's surface finish may result in loss of proper mechanical function of the device.
- Never attempt to reuse. Once the Eclipse Soft Tissue Anchor has been removed from the packaging, the device should be either discarded or sterilized. Never attempt to reuse the implant, even though it may appear undamaged.
- The surgeon should make the final decision regarding implant removal.
- In the absence of a bursa or joint, removal of the implant in elderly or debilitated patients is not recommended. Extreme care must be taken when removing the device.
- Use only Eclipse Soft Tissue Anchor components. Eclipse devices manufactured by MedShape must not be used in conjunction with other devices manufactured by any other manufacturer, as component parts may not be compatible.

6. POTENTIAL ADVERSE EFFECTS

- Potential adverse events resulting from use of the Eclipse Soft Tissue Anchor include, but are not limited to, the following:
- Loosening, cracking or fracture of the implant components.
- Loss of fixation in bone.
- Deep or superficial infection.
- Irrational injury of soft tissues.
- Sensitivity, allergies, or other reaction to the device material.
- Tissue reactions including macrophage and foreign body reactions adjacent to implants.
- Pain, discomfort, or abnormal sensations due to presence of the implant.
- Hematoma or thrombosis.

Adverse effects may necessitate re-operation, revision or removal surgery. Implant removal should be followed by adequate postoperative management.

7. PATIENT SELECTION INFORMATION

- The surgeon is responsible for patient selection.
- The surgeon is responsible for understanding the appropriate indications and contraindications associated with the device and selecting the surgical procedures and determinants to be best for each individual patient. Each surgeon must evaluate the appropriateness of the device and the procedure used to implant the device based on his/her own training and experience.
- The surgeon must determine if the device is appropriate for patients having any of the following physical or emotional conditions:
 - Drug and/or alcohol and/or smoke addiction and/or abuse.
 - Infectious disease.
 - Malignancy.
 - Local bone tumors.
 - Systemic or metabolic disorders.
 - Compromised wound healing.
 - Obesity.
 - Demonstrated psychological instability, inappropriate motivation or attitude.
 - Unwillingness to accept the possibility of multiple surgeries for revision or replacement.
 - Lacks an understanding that their preoperative capacity may not be fully recovered even after successful implantation.

8. PATIENT COUNSELING INFORMATION

- It is the responsibility of the surgeon to provide the patient with appropriate information prior to surgery. The surgeon should discuss with the patient all potential risks versus benefits of the procedure and the need for postoperative care and the expectations for improvement of his/her condition postoperatively. The patient should not have unrealistic expectations regarding the results that the surgery and implant may provide. In order to make an informed decision, the patient should clearly understand all applicable warnings, precautions, possible intraoperative and postoperative complications, and possible adverse effects associated with the surgical procedure and implantation of the device.
- The patient should be provided with detailed postoperative instructions regarding postoperative care, and the limitations of the device. Postoperative care and physical therapy should be structured to prevent excessive loading of the operative extremity until sufficient healing has occurred. The patient should be advised that noncompliance with postoperative instructions could lead to loss of fixation or device failure requiring revision surgery to remove the device. The patient should be encouraged to report to his/her surgeon any unusual changes to the operated extremity. If evidence suggests fixation failure, breakage, or migration of the implant, an intensified schedule of check-ups is advised and new warnings and instructions to the patient may be necessary to further restrict activities.
- The patient should be encouraged to receive prompt medical attention for any infection that may occur, either at the surgery site or elsewhere in the body.

9. PREOPERATIVE PLANNING INFORMATION

- Careful preoperative planning must be conducted.
- Never attempt a surgical procedure with defective, damaged, or otherwise compromised instruments or implants. Inspect all components preoperatively to ensure that the device components and instruments are appropriate for use.
- It is the physician's responsibility to determine the correct size of Eclipse Soft Tissue Anchor to be implanted. The physician should always have an inventory of sterile implants on hand at the time of surgery to ensure availability of the optimum size for the patient. If any of the devices are damaged during attempted placement, additional devices of the same size should be available.
- Alternative fixation methods should be available for use in the event that the Eclipse Soft Tissue Anchor cannot be successfully implanted.

Handling of the Eclipse Soft Tissue Anchor and ancillary surgical instruments must be performed in accordance with aseptic handling practices to maintain sterility following sterilization by the manufacturer (Device) or by the healthcare facility (Instruments).

10. GENERAL DIRECTIONS FOR USE

- Carefully determine the diameter of the tendon, ligament, or graft to be fixedated using the MedShape graft sizing block provided in the tenodesis procedure pack.
- Prep an appropriately sized bone tunnel per standard procedure.
- Push tendon into tunnel using tendon fork. Optional guide pin is provided in procedure pack to ensure the tendon remains secure in the tunnel during device insertion. OR
 - 3a. Whip stitch the end of the graft using a robust suture (#2 non-absorbable suture). Ensure that each suture limb is at least eight (8) inches long after suturing.
 - 3b. Pass the Suture Lasso through the implant body and pass one suture limb through the Suture Lasso loop. For the small diameter implants, there is already available suture lasso to encompass the optional pre-loaded suture loop may be used to facilitate the passage of the whipped stitched suture strands through the bone tunnel.
 - 3c. Pull the suture limb tight so that the tendon is held securely at the tip of the Sheath. Secure both suture limbs in the suture cleat on the Gun.
- 4. Disengage the Knob Lock to rotate the implant body for desired orientation. Reengage the Lock.
- 5. Place the Sheath assembly into the bone tunnel and carefully insert, ensuring that the compressed portion of the Sheath is adjacent with the tendon.
- 6. Ensure that the Sheath is fully inserted in the desired orientation until the back of the device is flush with the cortical layer. The Strike Pin located on the back of the Gun may be gently tapped to facilitate insertion.
- 7. Squeeze the Trigger and disengage the Trigger Lock to release Trigger in preparation for deployment.
- 8. Squeeze the Trigger to advance the Bullet along the Retention Rod.
- 9. Once the Bullet begins to engage the Sheath, ensure that the Bullet is centralized within the Sheath. Continue to advance the Bullet by squeezing the Trigger until the Bullet is fully seated within the Sheath cannulae.
- 10. Fully release the Trigger.
- 11. Disengage Knob Lock and remove the Retention Rod by rotating counter-clockwise until fully disconnected from Sheath.
- 12. Discard the Gun.
- 13. Confirm deployment of the Eclipse Soft Tissue Anchor and adequate graft fixation.
- 14. Trim the excess suture if present. Complete surgical procedure as required.

11. HOW SUPPLIED
The Eclipse Soft Tissue Anchor is provided sterile for single use only. Carefully inspect sterile packaging for damage prior to use. If the sterile packaging is found to be damaged or open, do not use the device or attempt to resterilize. Call your MedShape sales representative for a replacement.

12. STORAGE
Store the Eclipse Soft Tissue Anchors in a dry place at room temperature (20°C to 25°C).

13. WARRANTY INFORMATION
Limited Liability:
Each Eclipse Soft Tissue Anchor is guaranteed for materials, function, and workmanship for a single patient use. MedShape shall not be liable, expressly or implied, for any damage which might arise or be caused, whether by the customer or by any of the users of the product, as a result of:

- a) Misuse, mishandling, and/or improper organization.
- b) Repairs or modifications performed other than by MedShape or a MedShape authorized repair facility.
- c) Use in any manner or medical procedure other than those for which it is designed; and any special, indirect, and/or consequential damages of any kind however caused arising from the sale or use of the product.

THE WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON MEDSHAPE'S PART

Return Conditions:
In the event the device must be returned for any reason, return the product in the original packaging. Contact Customer Service or an authorized MedShape representative to receive a return authorization number prior to returning shipment.

14. SYMBOLS

	Consult Instructions For Use. Read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury.
	The use of this product is intended to be limited to persons trained in the procedure and knowledgeable of the inherent risks. The patient should be fully informed about the need to limit activity during healing.
	Catalog Number.
	Lot Number.
	Quantity.
	Expiration Date.
	Do Not Reuse.
	Manufacturer.
	Do Not Use If Package Is Damaged.
	Sterilized Using Ethylene Oxide (ETO).
	MR Safe: Device poses no known hazards.
	CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.

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