


BEFORE USING THE DEVICE, PLEASE READ THE **INSTRUCTIONS FOR USE** COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:

The user should be able to: Read and understand the directions, warnings, and cautions.

BRACE REMOVAL:

- Detach quick release buckles on waist strap (i) and front of sling (iii), remove forearm strap (ii) and thumb strap (iv). (Fig A) Gently remove shoulder immobilizer.




RE-APPLICATION:

For ease of application, place the brace on a couch or bed with the shoulder strap/ring and waist strap fully extended and away from the brace. The wider, open end of the shoulder ring should be facing up. Make sure the sling is positioned so you can place the affected arm/shoulder into the arm sling envelope easily when seated on bed/couch. (Fig B)


Step 1. Apply Arm Sling

Place forearm into arm sling envelope as far back as possible. Secure thumb strap (iv) between thumb and forefinger at front of sling. Secure forearm strap to sling and cushion (ii-Fig. C)



Step 2. Apply Shoulder Ring

Using your unaffected arm, reach back and place your unaffected arm through the shoulder ring. Make sure the hook and loop closure of shoulder ring is facing up and away from body. (Fig. D) Secure the shoulder strap buckles at front of arm sling. (Fig. E)



Step 3. Apply Waist Belt

While still seated, make sure cushion on injured side is at waist/hip level and elbow is fixed at 90°. Narrow end of cushion should be positioned to back of forearm/elbow. Bring waist belt around back and attach to quick-release buckle (i) (Fig. F1) at front of cushion. Adjust waist belt as necessary for proper fit. (Fig.F2)




NOTE: The Swivel D-Ring Strap Tabs on the front and back of Shoulder Ring should be positioned high enough on the ring (approx. 1" - 2" from top of shoulder) to maintain arm in proper position and for ideal comfort during wear. See G1 (front), G2 (back) (Fig. G)

CONTRAINDICATIONS: None.

WARNINGS AND CAUTIONS:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.



CLEANING & STORAGE INSTRUCTIONS:


Hand wash in water (86°F/30°C), using mild soap. Air dry.

NOTE: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

Do NOT iron. Do NOT tumble dry. Do NOT bleach.

MATERIAL CONTENTS:

Nylon, Polypropylene, Polyurethane Foam, Cotton, Laminate Foam/Tricot, Polyacetal, Polyester, Acetal.



WARRANTY:

DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

INTENDED FOR SINGLE PATIENT USE.

NOT MADE WITH NATURAL RUBBER LATEX

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569
5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A.
DJOGlobal.com