ProCare® Quick-Fit™ W.T.O.

Patient Application Instructions



BEFORE USING THE DEVICE, PLEASE READ THE **INSTRUCTIONS FOR USE** COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE / INDICATIONS:

Treatment of deQuervains syndrome; Scaphoid injuries; Gamekeeper's thumb; and Wrist and thumb sprains and strains.

CONTRAINDICATIONS:

None

APPLICATION INFORMATION:

- Release single pull strap and thumb straps. Position hand in brace with cotton stockinet liner on top of hand. Thumb should rest on the abduction stay. Make sure the Palmer stay is positioned comfortably below the palmer crease and aligned with forearm.
- Wrap single pull strap around wrist area and secure on brace. Strap should be snug but not constrict.
- Thumb straps should wrap in opposing directions and secure on brace. Straps should be snug but not constrict.

WARNINGS AND PRECAUTIONS:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately. Never wear support over an open wound.

CLEANING INSTRUCTIONS:

Hand wash in warm water using mild soap, rinse thoroughly. AIR DRY.

Note: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

WARRANTY:

DJO, LLC will repair or replace all or part of unit and its accessories for material and workmanship defect for a period of six months from the date of sale.

FOR SINGLE PATIENT USE ONLY.

NOT MADE WITH NATURAL RUBBER LATEX.

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.











T 800.336.6569 D 760.727.1280 F 800.936.6569 2900 Lake Vista Drive I Lewisville, TX 75067 I U.S.A. DJOGlobal.com