ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INDICATIONS: The Procare MiniTrax Walker is designed to provide support and immobilization to the lower leg, ankle and foot. It may be suitable for the treatment of stable foot fractures, ligamentous injuries such as Plantar Fasciitis, tendon injuries such as Achilles tendon repair and lateral ankle sprains. Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICATIONS: Do not use if you are allergic to any of the materials contained in this product. This device is contraindicated for unstable fractures of the lower leg, ankle and foot.

WARNINGS AND PRECAUTIONS:

- This device is to be used under the supervision of a healthcare professional. The determination of when to apply the brace and the frequency and duration of use should be strictly at the discretion of the treating physician. Always consult with your physician before making changes to the brace.
 - Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity.
 - Do not use this device on patient's incapable of communicating physical discomfort.
 - Check your skin frequently for skin damage and skin irritation.
 - Patients without sensation (i.e. post-op anaesthesia, neuropathies, etc.) should be monitored frequently for skin damage, skin irritation or wound management.
 - Use caution when walking on slippery or wet surfaces to avoid injury.
 - Do not use over open wounds.
 - Do not tighten straps. This may result in reduction of blood flow or sensation.
 - Do not modify the device or use the device other than intended.
 - Do not use this device if it is damaged and/or the packaging has been opened.
 - If pain, swelling, changes in sensation or other unusual reactions occur while using this product, consult your medical professional immediately.
- NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INFORMATION:

A cotton or cast sock will enhance comfort - NOT INCLUDED

1) PREPARE BRACE

- Unfasten straps and open liner.

2) APPLY BRACE

- While seated, place foot in brace (Fig. A).
- Position heel in back of brace. Wrap liner around foot and then ankle (Fig. B).
- Secure straps from bottom to top (Fig. C). Tighten until snug and comfortable.

3) REMOVE BRACE

- To remove the brace, unfasten the straps, open the liner and remove the boot.

CLEANING INSTRUCTIONS: The foam liner may be hand-washed in 86°F/30°C water with mild detergent and air-dried. Do not put in dryer or use other heat source to dry.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

FOR SINGLE PATIENT USE ONLY.

NOT MADE WITH NATURAL RUBBER LATEX.
RX Only.

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INDICACIONES: The Procare MiniTrax Walker is designed to provide support and immobilization to the lower leg, ankle and foot. It may be suitable for the treatment of stable foot fractures, ligamentous injuries such as Plantar Fasciitis, tendon injuries such as Achilles tendon repair and lateral ankle sprains. Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICACIONES: Do not use if you are allergic to any of the materials contained in this product. This device is contraindicated for unstable fractures of the lower leg, ankle and foot.

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CLEANING INSTRUCTIONS:

The foam liner may be hand-washed in 86°F/30°C water with mild detergent and air-dried. Do not put in dryer or use other heat source to dry.

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Solo con receta médica.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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