ENGLISH
BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: The REACTION knee brace is designed to help reduce pain by shock absorption and by returning kinetic energy to the knee. Indications for use are mild to moderate patellofemoral pain or maltracking and mild to moderate osteoarthritis of the knee.

CONTRAINDICATIONS: None

WARNINGS AND PRECAUTIONS: If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

APPLICATION INFORMATION:

1. Pull on sleeve until comfortable and patella is centered in the hour glass shaped opening. (Wearing sleeve is optional.)
- 2a. Lay web over top of knee joint, centering web by aligning opening with patella (figure 2b).

3. Thread tabs on one side of back pad through inside openings on web and affix hook portions to loop face of back pad. (Note: Spacer material side of back pad should be placed toward leg and loop side away from leg. Outside of web has the logo, inside of web is dimpled.)
4. Thread loose tabs of back pad through available openings on web and then affix hook tabs to loop surface on back pad.

- 5a. Once all tabs are affixed to assure proper tightness, lift and pull both of lower tabs until snug.
- 5b. Repeat for top pair of tabs -brace should be snug and comfortable, do not over tighten.

CLEANING INSTRUCTIONS: Both the plastic parts and pads can be hand wash with mild soap in lukewarm water. Rinse both thoroughly. AIR DRY separately.

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