INSTRUCTIONS FOR USE - AMF TT CONES

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GEBRAUCHSANWEISUNG - AMF TT CONES

MODE D'EMPLOI - AMF TT CONES

INSTRUCCIONES DE USO - AMF TT CONES

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INSTRUCTIONS FOR USE - AMF TT CONES

Before using a LimaCorporate product, the surgeon should carefully study the following recommendations, warnings and instructions, as well as the most updated product-specific information (e.g.: product literature, surgical technique).

1. PRODUCT INFORMATION

The AMF (Anatomic Metaphyseal Fit) TT Cones, including femoral and tibial cones, are knee prosthesis components intended to be used to: • reinforce the medullary cavity of the tibia and femur; • fill a proximal tibia or distal femur defect that may result from the removal of a primary knee system; • provide support of the tibial baseplate or femoral component by means of bone cement.

They are used in combination with the Physica cemented tibial plate and the Physica PS femoral component and Multigen Plus CCK and H Systems, which are manufactured by LimaCorporate.

The components included in this leaflet must not be used with components from other systems or components from other manufacturers.

AMF TT Cones are intended for uncemented fixation to the bone and are fixed to the femoral and tibial implants using bone cement.

1.1. MATERIALS

| Components | Material |
|--|----------|
| AMF TT Cones | Ti6Al4V |
| Material Standards | |
| Ti6Al4V (ISO 5832-3 - ASTM F1472) | |

Some patients may be hypersensitive or allergic to the implant materials; this should be appropriately considered by the Surgeon.

1.2. HANDLING AND STORAGE

All devices are provided sterile and should be stored at ambient temperature (indicative range 0-50° C / 32-122° F) in their protective closed packaging in controlled rooms, protected from exposure to light, heat and

sudden changes in temperature.

Once the package is opened check that both the model and size of the implant correspond exactly to the description printed on the labels. Avoid any contact between the implant and objects or substances that can alter the sterile condition or surface integrity. A careful visual examination of each implant is recommended before use in order to verify that the implant is not damaged. **Components removed from the package should not be used if they are dropped or suffer other accidental impacts. Devices should not be modified in any way.**

The device's code and lot number should be recorded in the patient's case history by using the labels included in the components packaging. The disposal of medical devices is to be performed by the hospitals in conformity with applicable laws.

Re-use of previously implanted devices must be absolutely avoided.

Risks associated with the reuse of single use devices are: • infection; • early or late failure of the device or device fixation; • lack of appropriate coupling between modular junctions (e.g. taper connections); • device wear and wear debris associated complications; • transmission of diseases (e.g. HIV, hepatitis); • immune system response / rejection.

The components forming original LimaCorporate systems must be rigorously coupled following the given indications. The use of single components, or components belonging to other systems, are subject to the approval by LimaCorporate. The manufacturer and the dealer are not liable for possible coupling incompatibility.

2. GENERAL INFORMATION ON INDICATIONS, CONTRAINDICATIONS AND RISK FACTORS

2.1. INDICATIONS

AMF TT Cones are intended for use in skeletally mature patients with bone defect or poor bone quality (osteoporotic bone) or in case of sclerotic bone that requires supplemental metaphyseal fixation in the clinical judgment of the surgeon. Tibial and femoral cones are intended for uncemented fixation to the bone and are fixed to the femoral and tibial implants using bone cement. Please also refer to the indications for use of total knee prosthesis (Physica system and Multigen Plus system) to be used in combination with AMF TT Cones.

2.2. CONTRAINDICATIONS

Absolute contraindications include: • acute or chronic, local or systemic infections; • minimal metaphyseal bone loss which not requires supplemental metaphyseal fixation in the clinical judgment of the surgeon; • insufficient / inadequate bone mass or quality which prevents a stable anchor of the implant in the clinical judgment of the surgeon; • progressive tumor diseases; • known incompatibility or allergy to the product materials; • septicemia; • persistent acute or chronic osteomyelitis; • open epiphyses (immature patient with active bone growth).

The relative contraindications are: • vascular or nerve diseases affecting the concerned limb; • metabolic disorders which may impair fixation and stability of the implant; • any concomitant disease and dependence that might affect the implanted prosthesis; • metal hypersensitivity to implant materials; • important osteoporosis,

haemophilic disease; • internistic problems with high risk for surgery; • skeletal immaturity.
Please also refer to the contraindications of total knee prosthesis (Physica system and Multigen Plus system) to be used in combination with AMF TT Cones.

2.3. RISK FACTORS

The following risk factors may result in poor results with this prosthesis: • overweight¹; • strenuous physical activities (active sports, heavy physical work); • fretting of modular junctions; • incorrect implant positioning (e.g. varus positioning); • wrong size of components; • improper fixation of the device; • insufficient bone to support the tibial or femoral components; • medical disabilities which can lead to an unnatural gait and loading of the knee joint; • muscle deficiencies; • multiple joint disabilities; • refusal to modify postoperative physical activities; • patient's history of infections or falls; • systemic diseases and metabolic disorders; • local or disseminated neoplastic diseases; • drug therapies that adversely affect bone quality, healing, or resistance to infection; • drug use or alcoholism; • marked osteoporosis or osteomalacia; • patient's resistance to disease generally weakened (HIV, tumour, infections); • severe deformity leading to impaired anchorage or improper positioning of implants; • use in combinations with products, prosthesis or instruments of another manufacturer; • errors of operative technique.

3. WARNINGS

3.1. PRE-OPERATIVE PLANNING

LimaCorporate products must be implanted only by surgeons familiar with the joint replacement procedures described in the specific surgical techniques.

ALLOWED/NOT ALLOWED COMBINATIONS

- Physica PS femoral component can be coupled with all sizes of the femoral cones;
- When using the bicondylar femur cones, the Physica PS femoral component must not be coupled with the optional PS pegs (6515.09.900);
- All sizes of tibial cones are compatible with Physica cemented tibial plate;
- When using Physica tibial stem, do not use the tibia cone 15;
- Multigen H tibial component cannot be coupled with central Tibia/Tibia cone size 18;
- Multigen H and CCK tibial component must not be coupled with Tibia cone size 15;
- Tibia/Central Tibia cone size 18 cannot be coupled with Multigen tem Dia. 20, 22, 24;
- Tibia/Central Tibia cone size 18 cannot be coupled with Multigen tibial module offset +3mm and +6 mm;
- Tibia/Central Tibia/ Peripheral Tibia cone size 21 cannot be coupled with Multigen stem Dia. 22, 24;
- Tibia/Central Tibia/ Peripheral Tibia cone size 21 cannot be coupled with Multigen tibialmodule offset + 3 mm and +6 mm;

¹ According to the definition of the World Health Organization (WHO), Body Mass Index (BMI) greater than or equal to 25 Kg/m²

- Tibia/Central Tibia/ Peripheral Tibia cone size 24 cannot be coupled with Multigen tibial module offset +3 mm, coupled with Multigen stem Dia.20, 22, 24;
- Tibia/Central Tibia/ Peripheral Tibia cone size 24 cannot be coupled with Multigen tibial module offset +6 mm;
- Tibia/Central Tibia/ Peripheral Tibia cone size 27 cannot be coupled with Multigen tibial module offset +3 mm, coupled with Multigen stem Dia. 22 and 24;
- Tibia/Central Tibia/ Peripheral Tibia cone size 27 cannot be coupled with Multigen tibial module offset +6 mm;
- Tibia/Central Tibia/ Peripheral Tibia cone size 30 cannot be coupled with Multigen tibial module offset +6 mm, coupled with Multigen stem Dia. 20, 22 and 24;
- Multigen H and CCK femoral component must not be coupled with Central Femur cone size 15;
- Central Femur cone size 18 cannot be coupled with Multigen femoral module R-L, coupled with Multigen stem Dia. 20, 22 and 24;
- Central Femur cone size 18 cannot be coupled with Multigen femoral module R+3mm/L-3mm or Multigen femoral module R-3mm/L+3mm, coupled with Multigen stem Dia. 20, 22 and 24
- Central Femur cone size 21 cannot be coupled with Multigen femoral module R-L, coupled with Multigen stem Dia.22 and 24;Central Femur cone size 21 cannot be coupled with Multigen femoral module R+3mm/L-3mm or Multigen femoral module R-3mm/L+3mm, coupled with Multigen stem Dia. 22 and 24;
- Bicondylar Femur cone size 18 cannot be coupled with Multigen femoral module R-L, coupled with Multigen stem Dia. 20, 22 and 24;
- Bicondylar Femur cone size 18 cannot be coupled with Multigen femoral module R+3mm/L-3mm or Multigen femoral module R-3mm/L+3mm, coupled with Multigen stem Dia 20, 22 and 24;
- Bicondylar Femur cone size 21 cannot be coupled with Multigen femoral module R-L, coupled with Multigen stem Dia. 22 and 24;
- Bicondylar Femur cone size 21 cannot be coupled with Multigen femoral module R+3mm/L-3mm or Multigen femoral module R-3mm/L+3mm, coupled with Multigen stem Dia , 22 and 24;
- Central Femur/ Bicondylar Femur cone size 24 and 27 must not be used when Multigen femoral module R+3mm/L-3mm or femoral module R-3mm/L+3mm are used.

Pre-operative planning, through radiographic templates in different formats, provides essential information regarding the type and size of components to be used and the correct combination of required devices based on the anatomy and specific conditions of each patient. Inadequate pre-operative planning can lead to improper selection of the implants and/or incorrect implant positioning.

Complications or failures of the total knee replacement are more likely to occur in heavy and highly active patients. The surgeon must perform a careful evaluation of the patient's clinical condition and level of physical activity before performing surgery.

LimaCorporate specialized technical staff is available to provide advice regarding pre-operative planning, the surgical technique, and product and instrumentation assistance both prior to and during surgery.

The patient must be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, that it has a finite expected implant life, and

may need to be replaced at some time in the future.

The possible impact of the factors mentioned in sections 2 and 3.4 should be considered preoperatively and the patient informed as to what steps he/she can take to reduce the possible effects of these factors.

Implants are single use devices; do not re-use implants that were previously implanted in another patient. Do not re-use an implant that has previously come into contact with the body fluid or tissue of another person.

Surgical instruments are subject to wear with normal usage. After extensive use or excessive loads, instruments are susceptible to fracture. Surgical instruments should be used only for their specific purpose. Before use, the functionality of surgical instruments should be checked since the use of damaged instruments may lead to early failure of the implants. Damaged instruments should be replaced before surgery.

3.2. INTRAOPERATIVE

The use of trial devices is recommended to check the correct site preparation, size and positioning of the implants to be used. It is recommended that additional implants are available during surgery for use in those cases requiring prostheses of different sizes or when the preoperatively selected prostheses cannot be used. The correct selection, as well as the correct seating/placement of the implant, is extremely important. Improper selection, positioning, alignment and fixation of the implant components may result in unusual stress conditions that may affect system performance and survival of the implant.

The surgeon must check that no cement penetrates between the cone and the bone during implantation; the cement shall be used only between the knee implant and the cone.

The components forming original LimaCorporate systems must be assembled according to the surgical technique and used only for the labeled indications. Use only instruments and prosthesis trials specifically designed for use with the implants being used. The use of instruments from other manufacturers or the use of instruments designed for use with other systems can lead to inappropriate preparation of the implant site, incorrect positioning, alignment, and fixation of the devices followed by loosening of the system, loss of functionality, reduction of the durability of the implant, and the need for further surgery.

Care must be taken when handling Trabecular Titanium devices; they should not get in contact with fabrics that can release particles inside the trabecular structure.

3.3. POST-OPERATIVE CARE

Adequate post-operative care should be provided by the surgeon or other suitably qualified medical staff. Regular postoperative x-ray follows up is recommended to detect possible changes in position or condition of the implant or surrounding tissues.

The surgeon should make the patient aware of the limitations of limb function after knee arthroplasty and that caution is needed, especially in the post-operative period. Excessive physical activity or trauma to the involved knee can lead to premature failure through loosening, fracture, or abnormal wear of the prosthetic implants. The patient should be cautioned by the surgeon to govern activities accordingly and advised that the implants

may fail due to excessive joint wear.

In particular the following precautions should be presented to the patient by the surgeon: • avoid repeated high weight lifting; • keep body weight under control, overweight conditions may adversely affect the outcomes of joint replacement; • avoid sudden peak loads (consequences of activities, such as running and skiing) or movements which can lead to sudden stops or twisting; • avoid positions that can increase the risk of dislocation.

Lack of appropriate post-operative rehabilitation instructions and care can negatively influence the outcome of the surgical procedure.

3.4. POSSIBLE ADVERSE EFFECTS

Adverse effects that can occur in knee arthroplasty include: • loosening of the prosthetic components, for example due to osteolysis; • prosthesis dislocation and instability; • breakage or damage to the prosthetic implant; • instability of the system because of inadequate soft tissue balancing; • dissociation due to incorrect coupling of the devices; • subsidence; • infection; • tissue reaction to the implant material or abrasion; • local hypersensitivity; • local pain; • periprosthetic fractures, including intra-operative fractures; • bone resorption; • temporary or permanent nerve or tendon damage; • flexion contracture; • limited range of motion; • lengthening or shortening of the leg; • additional surgery.

Some adverse effects can lead to death. General complications include venous thrombosis with/without pulmonary embolism, cardiovascular or pulmonary disturbances, hematomas, systemic allergic reactions, systemic pain, blood loss and temporary or permanent paralysis.

4. STERILITY

a. Implants

All implantable components are provided sterile with a Sterility Assurance Level (SAL) of 10^{-6} . AMF TT Cones are sterilized by radiation or EtO. Do not use any component from a package that has been previously opened or appears to be damaged. Do not use implants after the expiration date printed on the label.

b. Instruments

Instruments are supplied non sterile and must be cleaned, disinfected and sterilized before use according to appropriate validated methods (refer to the "Instrument Care, Cleaning, Disinfections and Sterilization" brochure for validated sterilization parameters; this brochure is available on request or downloadable from www.limacorporate.com in the Products section). Users should validate their specific cleaning, disinfection and sterilization processes and equipment.

5. MAGNETIC RESONANCE IMAGING (MRI)

AMF TT Cones have not been evaluated for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. Since these devices have not been tested, LimaCorporate cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy. The risks associated with a passive implant in an MR environment have been evaluated and are known to include heating, migration and image artifacts at or near the implant side.

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Limacorporate S.p.A.