Fresh Stored Osteochondral Allografts

DESCRIPTION

RTI Biologics, Inc. (RTI) makes available fresh stored osteochondral (OC) allografts to implanting physicians. These fresh stored osteochondral allografts have a limited storage time in order to maintain viable chondrocytes. Fresh OC grafts are processed from donated human tissue and released using expedited release criteria to maintain graft utility. Final culture results and donor autopsy report (if an autopsy was performed) may not be available prior to release of the tissue and therefore may not be reviewed during the initial assessment of fresh stored osteochondral allografts. RTI believes that the implanting physician, following appropriate discussion with their patient, is the most qualified individual to make decisions regarding treatment options, including the use of fresh osteochondral allografts.

DONOR SCREENING AND TESTING (SUMMARY OF RECORDS)

LOT This symbol on the outer label indicates the unique number assigned to the tissue donor.

The donated human tissue utilized for this implant was recovered from a donor screened for risk factors associated with infectious diseases and medical conditions that rule out donation. The donor's blood was tested for relevant communicable diseases in a laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) or equivalent and registered with the U.S. Food and Drug Administration (FDA) for donor testing. The following test criteria were met for this donor:

| REQUIRED INFECTIOUS DISEASE TESTING | |
|--|-------------------------|
| BLOOD TEST | ACCEPTABLE RESULTS |
| HIV-1/ HIV-2 Antibody | Negative / Non-Reactive |
| Hepatitis C Virus Antibody | Negative / Non-Reactive |
| Hepatitis B Surface Antigen | Negative / Non-Reactive |
| Hepatitis B Core Antibody (Total) | Negative / Non-Reactive |
| Syphilis | Negative / Non-Reactive |
| Human T-Cell Lymphotropic Virus I/ II Antibody | Negative / Non-Reactive |
| HIV-1 / HCV NAT-TMA | Negative / Non-Reactive |

A licensed physician for RTI Biologics Inc. determined that the donor met eligibility requirements. The physician utilized available relevant information which may have included, but was not limited to: family/next-of-kin interview, medical/hospital records, donor

physical assessment, infectious disease test results, radiology/pathology reports, death certificate and autopsy report (if performed).

WARNINGS AND PRECAUTIONS FOR IMPLANTING PHYSICIANS

Unlike other RTI allografts, final culture results and donor autopsy reports (if an autopsy was performed) may not be available prior to release of the tissue and therefore may not be reviewed during the initial assessment of fresh stored osteochondral allografts. Information may become available following release that suggests the presence of potentially transmissible disease processes. RTI will continue to collect and review information pertaining to the suitability of fresh stored osteochondral allografts as it becomes available. If it is determined upon review by the RTI Medical Director that this information would be relevant to the release status of the osteochondral allograft, RTI will immediately forward the pertinent information to the implanting physician at the location on file.

Surgeons should discuss these potential adverse events with their patients. Because donor screening methods are limited, certain diseases may not be detected. The following complications of tissue transplantation may occur:

- Transmission or causation of diseases of unknown etiology and characteristics;
- Transmission of known infectious agents including, but not limited to, HIV, Hepatitis B, Hepatitis C, syphilis and bacteria;
- Immune rejection of transplanted grafts; or,
- Loss of function and/or integrity of transplanted tissue due to resorption, fragmentation, and/or disintegration including, but not limited to, associated loss of continuity, displacement, bending and/or fracture.

All of the same medical/surgical conditions that may complicate any surgical procedure may occur during or following the transplantation of an allograft.

Neither RTI nor its affiliates expressly or by implication WARRANT any tissue selected, recovered, processed or distributed by itself, distributors, representatives or agents. Please promptly report possible adverse outcomes to RTI.

TISSUE PROCESSING REAGENTS

The graft may retain residual amounts of the following reagents used during tissue processing: ascorbic acid, primaxin (imipenem with cilastatin), polymyxin B, povidone-iodine, saline or vancomycin.

The graft is packaged in a solution of Ham's F-12 Medium with 10% Fetal Bovine Serum (non-U.S. Origin). Traces of these fluids may remain on the tissue.

PACKAGE INSPECTION CRITERIA

Upon arrival or usage of product check the packaging and verify the following:

- Package is not torn or damaged.
- Product received matches order.
- Package is not leaking fluid.
- Product has not reached expiration date.

If any of the above are found (i.e. package is damaged, passed its expiration date, or is leaking fluids) **DO NOT USE THE TISSUE. If there are any questions or concerns contact Customer Service immediately.**

STORAGE

Refrigerated

Upon arrival, immediately place the graft in a refrigerator maintained between 1 – 10°C.

DO NOT FREEZE TISSUE.

GRAFT USAGE

GENERAL INFORMATION

- Inner package contains fluid to keep tissue hydrated.
- Once opened, the graft must be used for the current procedure or discarded.
- Tissue is for single patient, single use only.
- Do not sterilize.
- Adhere to instructions provided for surgical techniques and with all instrumentation.
- This allograft and all packaging materials used by RTI are latex free.
- The outer plastic packaging acts as a non-sterile dust cover for shipping and storage.

DIRECTIONS FOR PREPARATION

1. Remove dust cover or box (*may require cutting open plastic dust cover*).
2. Remove Tissue Utilization Record (TUR) card and other items for later review or use.
3. Compare labels and verify information matches. *If information does not match, do not use graft.*
4. Use aseptic technique to remove the pouch containing the allograft.
5. Place inner package onto sterile field.
6. Inner package contains fluid to keep graft hydrated.
7. Cut open package with sterile scissors.
8. Remove allograft and place in sterile bowl.
9. Add sterile fluid to bowl. the fluid from the inner package or an isotonic solution may be used. DO NOT USE WATER.
10. Allograft must be submerged completely in solution and remain hydrated until use.
11. Destroy tissue that is not used in surgery.

TISSUE UTILIZATION RECORD (TUR card)

The TUR card is designed for tracking implanted grafts. All information provided is kept confidential and used for product tracking only. This card must be completed by the operating room staff or clinician after the surgery. For further instructions please see those listed on the TUR card. Return completed copy to RTI to address or fax provided on TUR card.

CUSTOMER RETURNS AND COMPLAINTS

If this tissue was not received or ordered from RTI directly please notify your distributor of any complaints or adverse events. For all tissue ordered directly from RTI's distribution center please use the contact information below for complaints and returns.

CUSTOMER SERVICE

RTI Biologics, Inc.
Alachua, Florida, 32615, U.S.A.
Tel: 386.418.8888 | 800.624.7238 | FAX: 386.462.5533
www.rtx.com.

CTO Number: RTI Biologics, Inc. 100053

Table 2: Definitions of label symbols:

| | | |
|---------------------------|------------------|----------------------------|
| | EXP or | |
| See instructions for use | Expiration date | Storage temperature limits |
| | QTY | |
| Do not reuse | Measurements | For prescription use only |
| SN | REF | LOT |
| Serial number (Tissue ID) | Catalogue number | Lot number (Donor ID) |

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RTI Biologics, Inc. : Numéro d'inscription CTO 100053