



Before using this brace please read the following instructions completely and carefully. Correct application and care is vital to the proper functioning of this brace.

INTENDED USE / INDICATIONS:

Treatment of wrist sprains, strains, and Carpal Tunnel Syndrome (CTS).

CONTRAINDICATIONS:

None

APPLICATION INSTRUCTIONS:

1. Loosen all straps and position hand in brace. Make sure the palmer stay is positioned comfortably below the palmer crease and aligned with forearm.
2. Position adjustable dorsal stay casing to ensure maximum support and proper fit. Stay casing may be moved proximal or distal for desired controlled range-of-motion.
3. Use single pull strap to secure body of brace. Straps should be snug but not constrict.
4. Secure strap between thumb and forefinger.

CLEANING INSTRUCTIONS:

Hand wash in warm water using mild soap, rinse thoroughly. AIR DRY. Note: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

WARRANTY:

DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

WARNINGS AND PRECAUTIONS:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

PATENT PENDING

FOR SINGLE PATIENT USE ONLY

NOT MADE WITH NATURAL RUBBER

For more information, please contact your DJO Global Territory Manager or Authorized Distributor



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOglobal.com

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Information related to various health, medical, and fitness conditions and their treatment is not meant to be a substitute for the advice provided by a physician or other medical professional.