# Foot Cryo/Cuff<sup>™</sup> Gravity Patient Information

Patient:	Doctor:	

## What is the Cryo/Cuff Gravity?

The Aircast Cryo/Cuff Gravity combines cold and focal compression to better manage pain and reduce swelling. Continuous cold and focal compression are delivered using a simple, gravity-fed delivery system.

## How is the Cryo/Cuff Gravity used?

#### STEP 1 - Prepare the cooler (Figure 1)

- 1. Connect the blue tube to cooler.
- 2. Add cold water to line inside cooler.
- 3. Fill with ice to the line inside the cooler.
- 4. Lay insulation disk on top of ice. Attach cooler lid snugly.
- 5. Allow five minutes with occasional shaking to chill water.

# STEP 2 - Always apply <u>EMPTY</u> Foot Cryo/Cuff (Figure 2) (The first time, have another person help)

- Use only Aircast Cryo/Cuffs or ArcticFlow Cuffs with Aircast Cryo/Cuff Gravity Unit.
- Carefully and completely read instructions supplied with the specific Foot Cryo/Cuff prior to application. Please follow the instructions below to continue cooler preparation.

#### STEP 3 - Fill and pressurize Cryo/Cuff (Figure 3)

- · Connect blue tube to Cryo/Cuff.
- · Open cooler air vent.
- Raise cooler above the Cryo/Cuff until the cuff is full NOTE: To <u>AVOID</u> excessive pressure during use, <u>DO NOT</u> raise the cooler higher than 15" (38 cm) above the Cryo/Cuff.
- Close cooler air vent.
- Cooler can be disconnected by pressing quick-disconnect.








Figure 3.

### To re-fill cooler:

Rechill Cryo/Cuff as needed.

IMPORTANT: <u>EMPTY</u> the Cryo/Cuff <u>COMPLETELY</u> while still attached to the cooler.

When Cryo/Cuff is empty then empty the cooler. Repeat steps 1-5 in **STEP 1** to prepare the cooler for refilling.

#### To clean cooler:

The cooler can be cleaned with disinfectant wipes. The tube assembly can be cleaned with common mild detergent and must be completely submerged and hand dried.

#### To troubleshoot:

- Make sure unit has ice and water filled to the line indicated inside the cooler.
- 2. Make sure to fill cooler and Cryo/Cuff as instructions indicate.
- 3. Check Cryo/Cuff connection; make sure to listen for a click when connecting the hose to the Cryo/Cuff.
- Make sure the Cryo/Cuff is completely <u>EMPTY</u> before refilling.

### **⚠** CAUTION

- DO NOT use an elastic wrap with the Cryo/Cuff Gravity.
- Dressing used under the Cryo/Cuff should be applied lightly.
- Reduce pressure with any sense of discomfort, numbness or tingling of the limb.
- EMPTY the Cryo/Cuff after each use, and prior to refilling.
- Patient skin condition checks must occur every 1-2 hours on a routine basis.



FOR PRODUCT SUPPORT CALL 1-888-405-3251

