

PROCARE®

DESIGNED
FOR KIDS



MINITRAX

BIG BENEFITS
SPECIFICALLY DESIGNED FOR YOUR
SMALLEST PATIENTS

Ideal support and comfort for your smallest patients with mild to severe foot and ankle injuries


MiniTrax™ is the only pediatric walking boot to feature kid-friendly design elements that make wearing a boot less intimidating. Beneath the fun exterior though, the MiniTrax is a serious boot with big benefits.



Easy to apply and adjust with a simple strapping and liner system.



Lightweight design includes anatomically shaped plastic uprights that help to facilitate a better fit.



Low-profile rocker sole of the MiniTrax features a durable, non-slip tread tough enough to withstand the test of time—and of a kid.

COLOR CODED FOR
FAST AND ACCURATE
SIZE SELECTION



LARGE

YELLOW
GIRAFFE



MEDIUM

BLUE
WHALE



SMALL

ORANGE
OCTOPUS



All of this adds up to big benefits specifically designed for your smallest patients.



MiniTrax Sizing Guide							MiniTrax Replacement Liners	
Size	Part Number	Age in Years	Foot Length	Leg Length	US Shoe Size	Euro Shoe Size		
Small	79-95510	1 - 2.5	Up to 6" (15.2 cm)	Up to 9" (23 cm)	Up to 8.5	Up to 25	Small	79-95610
Medium	79-95520	2.5 - 5	6" - 7" (15.2-17.8 cm)	Up to 11" (28 cm)	8.5 - 11.5	25 - 29	Medium	79-95620
Large	79-95530	5 - 8	7" - 8" (17.8- 20.3 cm)	Up to 14" (36 cm)	11.5 - 1.5 (Youth)	29 - 33	Large	79-95630



DJO, LLC | A DJO Global Company
T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOGlobal.com/procare

Together in Motion™



DJO Global proudly sponsors:

