

REACTION WEB[™] KNEE BRACE

ACTION REACTION



WHAT'S YOUR REACTION?



Now available in vibrant colors! For a limited time only.

Get fast relief with a unique, elastomeric web that actually pulls pain away from your knee.

The REACTION WEB[™] is a responsive, webbed approach to knee pain that gives you a distinct alternative to the basic knee sleeve. Developed by a worldrenown orthopedic surgeon, Dr. William Nordt, REACTION WEB[™] can help decrease knee pain and increase function to get you back to your active lifestyle.

Indications: helps absorb shock and relieve knee pain caused by Chondramalacia, Patellar Tendonitis and Tendonosis, general Patellofemoral tracking issues, and mild OA.

Ordering Information: when ordering please choose the 'brace number + color'. Example: 11-0215-2-09000 would be a XS/S Reaction Knee Brace in Pink.

Unique and Different

- NEW! REACTION WEB in Vibrant Colors Available for a limited time.
- **NEW!** Better Fit and Comfort The new back pad and sleeve fit and conform better to the knee for a softer more breathable feel.
- Progressive Pain Relief The elastomeric web design helps to reduce pain by dispersing energy across the knee. The web acts like a spring to absorb shock and shifts the peak loads away from the painful area of the knee.

PDAC-assigned code of L1810.



Measurements listed below are the circumference taken around the thigh, $6^{\prime\prime}$ above mid-patella.

Reaction Knee Brace		Size	Measur	Measurements		
11-0215-2-XXXXX		XS/S	13"-18.5"	13"-18.5" (33-47 cm)		
11-0215-3-XXXXX		M/L	18.5"-23.	18.5"-23.5" (47-60 cm)		
11-0215-4-XXXXX		XL/XXL	23.5"-29	23.5"-29.5" (60-75 cm)		
11-0215-5-XXXXX		XXXL	29.5"-31"	29.5"-31" (75-79 cm)		
Replacement Under Sleeve						
11-0541-2		XS/S	13" -18.5'	13" -18.5" (33-47 cm)		
11-0541-3		Μ	18.5"-21"	18.5"-21" (47-53.25 cm)		
11-0541-4		L	21"-23.5"	21"-23.5" (53.25-59.5 cm)		
11-0541-5		XL	23.5"-26	23.5"-26.5" (59.5-67.25 cm)		
11-0541-6		XXL	26.5"-29	26.5"-29.5" (67.25-75 cm)		
11-0541-7		XXXL	29.5"-32	29.5"-32" (75-81.25 cm)		
Color Information = XXXXX						
Orange	14000	Red	01000	Green	11000	
Pink	09000	Yellow	08000	Blue	02000	
Grey	n/a (If color is not specified web will be gray)					



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569 1430 Decision Street I Vista, CA 92081-8553 I U.S.A.



00-2445 Rev A