



DONJOY[®]
VERSAROM[®] HIP



PRECISE, INDIVIDUALIZED PROTECTION AND SUPPORT.

The low-profile, [DonJoy® VersaROM® Hip](#) helps deliver effective, functional ROM restriction in 15° increments. Easy-to-use abduction/adduction push buttons and malleable hip and thigh shells help provide an individualized fit for increased patient comfort and compliance. The adjustable compression belt with easy-grip pull tabs anchors the hip joint, and the padded thigh cuff increases optimal alignment. With a universal fit, one brace fits right or left hip up to a 54" waist circumference.



Integrated belt compression with cinching straps:
Pull tabs anchor the hip joint and thigh cuff to promote optimal alignment

Malleable hip and thigh frames:
Conforms to patient's hips for optimal, individualized fit

Precise ROM Adjustments:
Range of Motion allows controlled hip flexion and extension adjustments in 15° increments

Telescoping Height Adjustment:
Helps to accommodate a wide range of heights and body shapes

Laterally placed abduction paddle with three-point pressure:
Helps to maintain the head of the femur in the acetabulum

Customized Fit:
Two optional push-button hinges adjust abduction and adduction for an individualized fit

Lightweight, slim-profile design:
Combined with the softgoods, helps to maximize patient comfort and compliance

Indications: DonJoy® VersaROM® Hip is indicated for use with post-op patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals.

Additional indications include:

- Hip instability
- Joint revisions or reconstruction
- Stable hip fractures
- ROM control
- After reposition of luxated hip joints
- Muscular imbalance
- Arthroscopic hip procedures
- Labral repair, with or without gluteus medius repair
- Post-op proximal hamstring repair
- Patients at risk to dislocate

Ordering information

Part number	Description
11-3333	DonJoy® VersaROM® Hip
11-3334	VersaROM® Hip & Cold Therapy Combo Kit



T 800.336.6569 | D 760.727.1280 | F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJGlobal.com/donjoy

Copyright © 2020 by DJO, LLC
MKT00-9388 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.