



SAFE. PRECISE. SMART.  
Cooling, Heating & Contrast Therapy



SAFE.  
PRECISE.  
SMART.

### MULTI-MODALITY TREATMENT IN ONE CONVENIENT, DIGITAL DEVICE

Kelvi™ provides cryotherapy, thermotherapy, and contrast therapy 100% digitally without any liquids—no ice, no water, no solutions, no mess of any kind.

KelviTEC™ technology measures the temperature of the skin, not the temperature of the device. By utilizing a continuous biofeedback loop to make continuous adjustments, Kelvi provides a safe, exact, measurable, and repeatable treatment—everytime.

### TREATMENT OPTIONS

Easy application in clinic, on the road, or at home.

#### CRYOTHERAPY

39°-59°F

#### THERMOTHERAPY

95°-108°F

#### CONTRAST THERAPY

Rapid comfortable transitions between cooling and heating


39°-108°F

### POST-OPERATIVE

- Cryotherapy has been shown to decrease pain and inflammation\*
- Physician/Clinician can customize programs to fit patient needs

### RECOVERY

- Extended treatment sessions for post activity needs
- Post-treatment rehab/cooling or contrast with any wrap
- Broad contrast therapy options without tub submersion
- Easy application in clinic, on the road, or at home



### CLINICIAN APP

Prescribe specific temperature and time for treatment protocols

Create and save customized protocols that can easily be pushed to the patient

### ON CONSOLE:

- Single temperature
  - Choose cold or hot
  - Select from 4 preset temperatures
- Contrast treatment
- Patient can adjust pressure to fit throughout treatment as needed



Rapid (~6sec) tile temperature change



Wi-Fi and Bluetooth compatible for easy software updates and App connectivity





NOW DISTRIBUTED BY DJO®

*The Kelvi™ PRO<sup>e</sup> is the smaller, portable unit in the Kelvi Pro series. With a simple, interactive display, and expanded features, the evolution of Kelvi technology is bringing more efficient Cooling, Heating, and Contrast therapy to home, clinic, and sports medicine facilities.*

DESCRIPTION	PART NUMBER
Kelvi™ PRO <sup>e</sup> Console	KLV10003
<b>ACCESSORIES</b>	
Kelvi PRO <sup>e</sup> Power Supply	KLV40004
Anti-Microbial Drops	KLV40002
Utility Strap	KLV50005
<b>SMART WRAPS</b>	
Back PRO (includes sleeve)	KLV20002
Knee PRO (includes sleeve)	KLV20003
Ankle PRO (includes sleeve)	KLV20004
Shoulder PRO (includes sleeve)	KLV20001
<b>SMART WRAP SLEEVES</b>	
Back PRO Sleeve	KLV30002
Knee PRO Sleeve	KLV30003
Ankle PRO Sleeve	KLV30004
Shoulder PRO Sleeve	KLV30001
<b>CASES</b>	
Smart Wrap Case Large	KLV50003
Smart Wrap Case Small	KLV50004
Kelvi PRO <sup>e</sup> Hard Case	KLV50007
Kelvi PRO <sup>e</sup> Soft Case	KLV50008
Smart Wrap Case Carrier	KLV50006

## ELEVATE YOUR TREATMENT WITH EMERGING TECHNOLOGY

Kelvi™ PRO<sup>e</sup>

Kelvi™ PRO



BACK PRO



KNEE PRO



ANKLE PRO



SHOULDER PRO

**enovis**

T 800.336.6569 F 800.457.4221

DJO, LLC  
5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A.  
enovis.com/kelvi

Copyright © 2022 DJO, LLC  
MKT00-12674 Rev A

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

kelvipro
 kelvipro
 kelvi\_pro

Kelvi™ is a trademark of Kelvi – Hypothermia Devices, Inc

**WARNING:** This device can be cold or hot enough to cause serious injury. Serious adverse reactions & safety hazards may occur when using this device. Read and understand all warnings and Instructions for Use, including contraindications, before using this device. Additional warnings appear in the Instructions for Use. Patients can only use this device with a prescription from a healthcare professional. The prescription must state how long and how often the device should be used and the length of breaks between uses.