



PROCARE[®]

MaxTrax™ 2.0

Engineered for customized healing



The road to healing has many curves. Navigate well with the new MaxTrax™ 2.0, built on years of expertise and valuable feedback from clinicians and patients alike. Our MaxTrax 2.0 Walking Boot is engineered with customizable malleable double uprights, a wide stable footbed, cross-strapping, pneumatic and non-pneumatic liner—and adds new innovations allowing for greatly enhanced protection, comfort and customization.

PERSONALIZED FIT.
TRUSTED SUPPORT.
ENHANCED COMFORT.

CUSTOMIZABLE MALLEABLE UPRIGHTS: hybrid plastic and metal construction accommodate post-trauma requisites

CROSS STRAPPING VERSATILITY: for customized fit

ACCOMMODATES A WIDER RANGE OF LEG AND FOOT SIZES with improved boot length and liner circumference

ADJUSTABLE DUAL AIR-CELL SYSTEM to maintain limb contact for improved stability and accommodate changes in edema

IMPROVED HEEL STABILITY: Achilles strap helps to secure heel during wear and may be removed as needed



CONVENIENT SIDE-MOUNTED PUMP and release valve on pneumatic styles

PROMOTES NATURAL GAIT AND SHOCK ABSORPTION with optimized, low-profile footbed and non-slip rocker sole



ADJUSTABLE STRAPS allow for repeatable application and flexibility to move around bandages or hot spots

CUSTOMIZABLE FIT AND IMPROVED COMFORT with unique molded foam insole. Raised rounded counter heel is fully trimmable to accommodate surgical sites.

MaxTrax™ 2.0 walking boots are available in pneumatic and non-pneumatic, tall and ankle heights.



MaxTrax™ 2.0 Air, Tall



MaxTrax™ 2.0 Air, Ankle



MaxTrax™ 2.0 Standard, Tall



MaxTrax™ 2.0 Standard, Ankle

MAXTRAX 2.0 AIR WALKER

TALL AIR	TALL AIR REPLACEMENT LINER	ANKLE AIR	ANKLE AIR REPLACEMENT LINER	MENS SHOE	WOMENS SHOE	KIDS SHOE	SIZE
79-95722	25-2112-1	79-95732	25-2113-1	<3.5	<4.5	<5	XS
79-95723	25-2112-2	79-95733	25-2113-2	3.5-7.5	4.5-8.5		S
79-95725	25-2112-3	79-95735	25-2113-3	7.5-10	8.5-11		M
79-95727	25-2112-4	79-95737	25-2113-4	10-14	11-15		L
79-95728	25-2112-5	79-95738	25-2113-5	14-18	15-19		XL

MAXTRAX 2.0 WALKER

TALL	TALL REPLACEMENT LINER	ANKLE	ANKLE REPLACEMENT LINER	MENS SHOE	WOMENS SHOE	KIDS SHOE	SIZE
79-95702	25-2110-1	79-95712	25-2111-1	<3.5	<4.5	<5	XS
79-95703	25-2110-2	79-95713	25-2111-2	3.5-7.5	4.5-8.5		S
79-95705	25-2110-3	79-95715	25-2111-3	7.5-10	8.5-11		M
79-95707	25-2110-4	79-95717	25-2111-4	10-14	11-15		L
79-95708	25-2110-5	79-95718	25-2111-5	14-18	15-19		XL

enovis

T 800.336.6569 F 800.457.4221

DJO, LLC
5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A.
enovis.com/procare

Copyright © 2022 DJO, LLC
MKT00-10583 Rev A

Indications: Support and management of edema for acute ankle sprains, soft tissue injuries of the lower leg, stress fractures of the lower leg, stable fractures of the foot and ankle, stable delayed union or non-union fractures of the distal tibia and fibula, Achilles tendon repairs, bunionectomies, metatarsal fracture, forefoot, midfoot and hindfoot injury.

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.