

**PROCARE**<sup>®</sup>

WHY **UNIVERSAL** IS A  
PERFECT FIT





## WHY UNIVERSAL IS A PERFECT FIT

*New universal designs deliver a custom fit, patient comfort and functional treatment for trauma or injuries to the wrist, hand and thumb.*



### QUICK-FIT W.T.O.

#### *Functional Treatment For Improved Outcomes*

A variety of wrist and thumb pathologies including DeQuervain's syndrome, Gamekeeper's thumb, and Scaphoid injuries may be addressed with this universal configuration that offers ease of application and contoured fit.

### QUICK-FIT WRIST AND FOREARM

#### *Universal Immobilization Post Injury*

A 10" length, adjustable dorsal stay pod and singular closure system provide an optimal immobilization solution for post-surgical or post-fracture casting of the wrist and forearm.

### QUICK FIT® WRIST II

#### *Versatility To Fit The Anatomy*

The simple, single-pull lace closure design of the Quick-Fit Wrist II and the adjustable, malleable dorsal stay pod provide variability and a contoured, custom fit to accommodate the changing wrist anatomy during rehabilitation, post-trauma or injury.

#### *Lace Closure System*

Simple, easy application is what you want for your patients to help improve compliance, outcomes and rehabilitation time. That's why universal is the perfect fit.



### UNIVERSAL THUMB-O-PRENE™

#### *Enhanced Design For Maximum Patient Comfort*

The newly enhanced Universal Thumb-O-Prene offers a wrap-around bilateral design with a low-profile, comfortable, custom fit that is both functional and easy to apply. Unique options allow for a variety of pathologies of the thumb to be addressed including stabilization of the CMC joint, Gamekeeper's thumb, and tendonitis.

### ProCare® THUMB SPLINT

#### *Protection With A Customized Fit*

The new ProCare Thumb splint provides complete stabilization of the thumb without restricting the use of the wrist – all in a universal support. Two stays, rigid and spiral provide needed immobilization for various thumb pathologies. Step-down rehab is an easy option.



## UNIVERSAL IS A PERFECT FIT FOR YOUR PRACTICE AND YOUR PATIENTS

### ORDERING INFORMATION

#### Quick-Fit® Wrist II

Patent pending, single-pull closure provides versatility to accommodate various patient sizes and rehabilitation phases. A soft, breathable liner helps maximize patient comfort and the adjustable, malleable dorsal stay pod may be used to customize immobilization to accommodate protocol. The unique contoured palmar stay offers proper anatomical fit and full finger dexterity. Patent pending.

Indications: Wrist sprains, strains, Carpal Tunnel Syndrome, post-wrist and hand trauma and rehabilitation.

Right	Left	Wrist Circumference	Length	Size
79-87560	79-87570	≤ 10"	8"	U
79-87561	79-87571	10" +	8"	U / XL

#### Quick-Fit Wrist & Forearm

Comfortable stabilization is provided with cool, soft construction and an easy single-pull lace closure for convenient application. An adjustable dorsal stay pod can be moved proximal or distal for desired controlled range-of-motion and the secondary dorsal stay may be removed to accommodate rehabilitation needs.

Indications: Post surgical immobilization, post fracture casting, wrist / forearm sprains and strains.

Right	Left	Wrist Circumference	Length	Size
79-87500	79-87510	≤ 13"	10"	U
79-87501	79-87511	13" +	10"	U / XL

#### ProCare® Thumb Splint

Unique perforated neoprene design with dorsal, radial and palmar stay pockets for a custom fit and maximum stabilization of the thumb. Each pocket contains a malleable aluminum stay and spiral stay that may be removed or shaped for desired support and custom fit.

Indications: Gamekeeper's/Skier's thumb, Basal Joint Arthritis and thumb tendonitis.

	Size
79-82710	U

#### Quick-Fit W.T.O.

Universal wrist, thumb orthosis with a single pull lace-closure design allows for ease of application and a secure, customizable fit delivering ease of use and functional treatment. Lightweight foam laminate maximizes patient comfort and the pre-contoured aluminum palmar and thumb stay may be adjusted for proper angulation to address wrist and thumb pathologies. Patent pending.

Indications: DeQuervain's Syndrome, Gamekeeper's thumb, scaphoid injuries, sprains and strains.

Right	Left	Wrist Circumference	Size
79-87480	79-87490	≤ 13"	U
79-87581	79-87591	13" +	U / XL

#### Universal Thumb-O-Prene™

Neoprene, bi-lateral wrap-around design offers two products in one with a low profile, custom fit. Two stays, rigid and spiral provide needed immobilization for various thumb pathologies or may be removed for step-down rehabilitation. The addition of unique, movable elastic strap provides compression to the base of the thumb to help stabilize the CMC joint.

Indications: Instability of the CMC (Carpometacarpal) joint, Basal Joint Arthritis and thumb tendonitis.

	Size
79-82700	U
79-82701	U / XL



T 800.793.6065 D 760.727.1280 F 800.580.9625  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
djoglobal.com

*Together in Motion™*

DJO Global proudly sponsors:

