

Designed to meet the needs of demanding athletes.

The Stabilizing Pro and new Stabilizing Speed Pro ankle braces reduce the chance of ankle injury by 61%¹. Both feature a new contoured design made with rugged ballistic nylon, reinforced eyelets, 3" circumferential compression strap and removable articulated stays. The Stabilizing Pro braces give trainers confidence that the vulnerable ankles of their athletes are protected. In addition to all the benefits of the Stabilizing Pro, the new Stabilizing Speed Pro brace features a "speed-lacing" system that makes applying the brace exceptionally fast and easy.

The Stabilizing Pro and Stabilizing Speed Pro: 61% fewer ankle injuries. 100% DonJoy.







ORDERING INFORMATION

Stabilizing Pro Ankle	Color	Size
11-3234-X-06000	Black	XXS-XXXL
11-3234-X-15000	White	XXS-XXXL
Stabilizing Speed Pro Ankle	Color	Size
11-3235-X-06000	Black	XXS-XXXL

Sizes Available (Measurement is circumference, taken from the base of the heel around ankle, while wearing an athletic sock)

X = 0 = XXS9" - 10" (23-25 CM)	X = 4 = L13" - 14" (33-36 CM)
X = 1 = XS10" - 11" (25-28 CM)	X = 5 = XL14" - 15" (36-38 CM)
X = 2 = S11" - 12" (28-30 CM)	X = 6 = XXL15" - 16" (38-41 CM)
$X = 3 = M \dots 12'' - 13'' (30-33 CM)$	X = 7 = XXXL16" - 17" (41-43 CM)

^{*}With athletic sock on, measure circumference of ankle from base of heel around ankle. If in doubt, size down to ensure maximum support.



DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569 1430 Decision Street I Vista, CA 92081-8553 I U.S.A. DJOGlobal.com/donjoy

Together in Motion...

References:

1. McGuine T.A., Brooks A. AM J Sports Med.2011 "The Effect of Lace-up Ankle Braces on Injury Rates in High School Football Players"

DJO Global proudly sponsors:







