DONJOY®



The X-Act ROM™ Elbow



Protect and Preserve. Precisely.

The X-Act ROM™ Elbow brace is a versatile, easy to apply, comfortable telescoping brace, allowing for ideal fixed or protected range of motion to help improve patient outcomes.

The shape-to-fit formable aluminum cuffs and plush thermoformed pads contour this lightweight, open design to the patient's elbow, providing a precise fit following surgery or injury. Telefit technology and snap-to-fit buckles allow for ideal adjustment and comfort during rehabilitation.





Easy to adjust telefit sliders for quick length adjustment and customized, precise fit.



The X-Act Rom Elbow brace ships flat for quick and easy application. Anti-migration material is integrated into thermoformed pads.

Swivel & Snap

Custom Fit

All four cuffs telefit to size

allowing for a more precise

fit on any size patient without impingement

Make it fit and make it fast with swivel action buckles that easily quick clip with one hand for a contoured fit



Precise Protection

Single, lateral aluminum hinge design helps to maximize protection and minimizes brace weight. User friendly, preferred flexion ROM settings from 10° to 120° in 10° increments. Easy quick-lock allows for fixed/locked position from 0° to 90° in 15° increments.

Easy, shape-to-fit malleable bicep and forearm cuffs contour to the anatomy for a custom fit and minimal application time.



Snap-to-fit buckles for faster fit with no hassle.

Maximum Comfort

Plush ComfortForm™ padding helps to maximize comfort for long-term wear



Supination/Pronation accessory provides the wrist with static stretch.





Ordering Info

Part # 11-9121	X-Act ROM Elbow, Right
Part # 11-9122	X-Act ROM Elbow, Left

Accessories

Part # 25-6753	Shoulder Sling
Part # 25-6757	ComfortForm Replacement Pad Kit
Part # 25-6758	Elbow Strap Replacement Kit
Part # 25-6755	Brace Lock Kit
Part # 25-6756	Pronation/Supination Kit, Right
Part # 25-6759	Pronation/Supination Kit, Left

Call 800.336.6569 or visit djoglobal.com



DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569 1430 Decision Street | Vista, CA 92081-8553 | U.S.A. DJOglobal.com

Copyright © 2018 by DJO, LLC





