

The Kelvi™ PRO^e is the smaller, portable unit in the Kelvi Pro series. With a simple, interactive display, and expanded features, the evolution of Kelvi technology is bringing more efficient Cooling, Heating, and Contrast therapy to home, clinic, and sports medicine facilities.

Kelvi™ PRO^e

Features:

- · All digital Cooling (39-59F) and Heating (95-108F) therapy
- · Rapid contrast therapy in a convenient, portable, and hygienic modality
- Kelvi™ wraps contain digital thermometers to measure skin temperature and create a comfortable and safe treatment
- · Biofeedback loop ensures continuous precise treatment
- · No changing ice and/or water; No solutions; No mess
- · Wraps do not mold from condensation

Why Kelvi™?

- Safe: Digital thermometers embedded in the KelviTEC™ tiles, constantly measure the patient's skin temperature
- Precise: Select your exact temperature and your exact time based on the patient's individual needs
- Smart: WiFi and Bluetooth connectivity expand your treatment options in this advanced system





Kelvi™ PRO^e

Description	Part Number
Kelvi™ PRO ^e Console	KLV10003
Accessories	
Kelvi PRO ^e Power Supply	KLV40004
Anti-Microbial Drops	KLV40002
Utility Strap	KLV50005
Smart Wraps	
Back PRO (includes sleeve)	KLV20002
Knee PRO (includes sleeve)	KLV20003
Ankle PRO (includes sleeve)	KLV20004
Shoulder PRO (includes sleeve)	KLV20001
Smart Wrap Sleeves	
Back PRO Sleeve	KLV30002
Knee PRO Sleeve	KLV30003
Ankle PRO Sleeve	KLV30004
Shoulder PRO Sleeve	KLV30001
Cases	
Smart Wrap Case Large	KLV50003
Smart Wrap Case Small	KLV50004
Kelvi PRO ^e Hard Case	KLV50007
Kelvi PRO ^e Soft Case	KLV50008
Smart Wrap Case Carrier	KLV50006

Kelvi.com **②** kelvipro **f** kelvipro **У** kelvi_pro

enovis...

T 800.336.6569 F 800.457.4221

DJO, LLC 5919 Sea Otter Place, Suite 200 | Carlsbad, CA 92010 | U.S.A enovis.com/donjou

Copyright © 2022 DJO, LLI MKT00-12590 Rev A **WARNING:** This device can be cold or hot enough to cause serious injury. Serious adverse reactions & safety hazards may occur when using this device. Read and understand all warnings and instructions for Use, including contraindications, before using this device. Additional warnings appear in the Instructions for Use. Patients can only use this device with a prescription from a healthcare professional. The prescription must state how long and how often the device should be used and the length of breaks between uses.

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advict The contents of this document do not constitute medical, legal, or any other type of professional advice Rather, please consult your healthcare professional for information on the courses of treatment, if any which may be appropriate for you.