



**REACTION WEB® KNEE**

<b>SIZING</b>			
Thigh Measurements taken 6" above mid-patella			
	PART NUMBER	THIGH CIRCUMFERENCE	SIZE
	11-0215-2-XXXXX	13" - 18½" (33 - 47 cm)	XS/S
	11-0215-3-XXXXX	18½" - 23½" (47 - 60 cm)	M/L
	11-0215-4-XXXXX	23½" - 29½" (60 - 75 cm)	XL/XXL
	11-0215-5-XXXXX	29½" - 31" (75 - 79 cm)	XXXL

**REDDIE™ BRACE**

<b>SIZING</b>			
Measurements taken 6" above mid-patella			
	PART NUMBER	THIGH CIRCUMFERENCE	SIZE
	79-82392	13½" - 15½"	XS
	79-82393	15½" - 18"	S
	79-82395	18" - 20½"	M
	79-82397	20½" - 23"	L
	79-82398	23" - 25½"	XL
	79-82399	25½" - 28"	XXL
	79-82399-10	28" - 30½"	XXXL

For product assistance, please contact Product Support at 1.888.405.3251 or email [product.specialist@djoglobal.com](mailto:product.specialist@djoglobal.com)

Please read the official **Instruction for Use** prior to using this device. Adjustments can only be made by a licensed professional.

Individual results may vary. Neither DJO, LLC nor any of its subsidiaries dispense medical advice. The contents of this poster do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.



DJO, LLC | A DJO Global Company

T 800.336.6569 D 760.727.1280 F 800.936.6569  
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.

[DJOGlobal.com](http://DJOGlobal.com)